

Clinitrend[®] Product Data Sheet

Product Description: L320 Latex Free Forehead Moving Line Temperature Monitor

Standard Range: 84-106°F / 29-41°C

Construction:

Adhesive Tape Adhesive: Acrylic copolymer pressure sensitive adhesive, thickness 37 μ each side of carrier
Carrier: Non-woven tissue, thickness 33 μ
Release Liner: Fully bleached, thickness 178 μ
Operating temp: -1 to 82°C (continuous use)

Covering Film Polyester

Color Change Material Nontoxic mixture of liquid crystal esters

Accuracy $\pm 1.00^\circ\text{F}$ / 0.6°C

Shelf Life 12 months from invoice date when stored at room temperature & humidity
(i.e. 21°C (70°F) & 50% relative humidity)

Application Data: Peel temperature-indicating label from backing paper. Apply to forehead. Reversible analog line moves along the scale to show detailed temperature readings in °F and °C.

Function: Continuous display of temperature, as indicated by a bright line, of anesthetized patients. Core adjusted with 3.5°F / 2°C offset.

Use: Ideal for use in Pre-Op, Surgery, PACU