

Britevent®

Urine Specimen Temperature Verification Thermometer

A comprehensive system of thermometers based on proprietary formulations of liquid crystals that change color at a designated temperature.



Product Code	Description	°C or °F	Range	Packaging
4030	Disposable Adhesive Label	°F	88-102°F	Sheets of 62 Kiss Cut (Die cut top layer but not backing adhesive paper) OR Individual Split Back (Easy to Remove Backing) bags of 100
4030B	Disposable Adhesive Label	°C	32-39°C	Sheets of 62 Kiss Cut (Die cut top layer but not backing of self-adhesive paper) OR Individual Split Back (Easy to Remove Backing) bags of 100



Product Code	Description	°C or °F	Range	Packaging
4100	Disposable Adhesive Label	°F	90-100°F	Sheets of 62 Kiss Cut (Die cut top layer but not backing of self-adhesive paper) OR Individual Split Back (Easy to Remove Backing) bags of 100
4486	Disposable Adhesive Label	°C	32-37°C	Sheets of 62 Kiss Cut (Die cut top layer but not backing of self-adhesive paper) OR Individual Split Back (Easy to Remove Backing) bags of 100

Product Highlights

- **Continuous Display Thermometer for Urine Specimen Collection Cup**
- **Easy to Read Digital Event Strip with Adhesive Backing**
- **Instant Response**
- **Conforms to SAMHSA Guidelines**
- **Sheets of 62 Kiss Cut (Die cut to top layer but not the backing of self-adhesive paper) OR Individual Split Back (Easy to Remove Backing) bags of 100**
- **Temperature Range: 88° - 102°F or 32° - 39°C (4030 or 4030B) 90° - 100°F or 32° - 37°C (4100 or 4486)**
- **Direct Read or 2°F Offset**

To Read:

1. Apply to cup.
2. Read Green Color for Temperature.

