

Brite Dot®

Specimen Temperature Verification Thermometer

Roll to Roll Labels for High-Volume Automated Application



Actual Size

Specifications

Product Code	D4630 ADR3000
Description	Adhesive backed label on a clear polyester release liner
Values	90-100°F & 32-38°C w/2°F Offset
Label Size	1.5" x .3937"
Packaging	Rolls of 3,000



To Read:

1. Apply to Cup
2. Read Green Color for Temperature

Product Highlights

- High-volume/Low-cost thermometer
- °F & °C Temperature Scale
- Appropriate for automated applications
- Adhesive backed roll stocked label
- Storage and handling convenience
- Temperature accuracy $\pm 1^\circ\text{F}/0.5^\circ\text{C}$