

Briteline

Product Data Sheet

Product Description: Briteline Urine Temperature Indicators

Standard Range: 89 to 101°F or 31.5°C to 38.5°C

Construction:

| | |
|-----------------------|--|
| Adhesive Tape | Adhesive: Single-sided adhesive paper label (Fasson S-246) |
| Covering Film | Polyester |
| Color Change Material | Non-toxic mixture of liquid crystal esters |
| Accuracy | L215F $\pm 1^{\circ}\text{F}$ L215C $\pm 1^{\circ}\text{C}$ L216F $\pm 1^{\circ}\text{F}$ |
| Warranty Period | 12 months from invoice date when stored in optimal storage conditions (between 60 and 80°F and 40-60% relative humidity) |

Application Data: Apply temperature indicating label to outside of urine specimen cup. Reversible analog line moves along the scale to show detailed temperature reading in °F or °C.

Function: Display of temperature, as indicated by a bright line, of urine specimen cup. Label temperatures are adjusted with a 1°F / 0.6 °C offset.

Use: Ideal for use in clinics, hospitals, and doctors' offices.